

Please ask me if I would like to spell a message, be asked a yes/no question or say a quick phrase. For spelling, follow along as I point to a letter. If I cannot point myself, slowly point to and say A, E, I, O, U. When I indicate* my desired letter is in the row, point to each letter in that row until I indicate* my desired letter. If I want to say a phrase, point to and say 1-6 until I get to my desired phrase.


get to my desired pindse.							
					(5) ESSENCE LOW TECH		
I want to spell it	Ask me a Y/N question	Go to my phrases	YES	NO	MAYBE	1	
Α	В	C	D	space	start over	2	
E	F	G	н	guess	mistake	3	
	J	К	L	Μ	Ν	4	
0	Р	Q	R	S	Т	5	
U	V	W	X	Y	Z	6	other side

*Using an eye blink is often the most reliable way for a person to signal their partner; however, other movements may be used so long as the partner knows which response to watch for while communicating with the individual (i.e. looking up, smiling, or tapping a finger).

prentrom.com

Essence[®] for Accent Devices By Prentke Romich Company





prentrom.com

Essence[®] for Accent Devices By Prentke Romich Company